


Filename: Observing exercise 4 – Web photo of Moon with Vision.docx

1. Start with the file "[Observing exercise 1](#) - "Go To" Saturn and Mars.pdf"
2. Set the telescope in "Home Position" and "Go To" the Moon
3. Put web camera in "90 eyepiece holder" and connect camera to USB port (Article 4 – Quick Start Guide – Image transfer directly from web camera to Internet)
4. Create a folder on your hard drive (C:\Metoichi Observations 03062014)
5. Load Vision (The program selects web camera automatically)
6. Go to "Device Settings" in "Devices Windows", choose "Archive" in the menu to left and open the folder (Store locally in directory) C:\Metoichi Observations 03062014.
7. We can change the camera settings: Click the icon on menu bar at the bottom of the screen to the right:  (the program: Philips ToUcam). Go to the tab "Bildekontroller" : "Full auto" (Article 6, page 14)